INSTRUCTIONS FOR USE

The NGAL Test™ Calibrator Kit

 ST002CA

INDEX

English	Page 2
Czech	Page 3
Danish	Page 4
French	Page 5
German	Page 6
Greek	Page 7
Italian	Page 8
Polish	Page 9
Russian	Page 10
Spanish	Page 11
Swedish	Page 12

EN
CZ
DA
FR
DE
GR
IT
PL
RU
ES
SE

The NGAL Test™

For your clinical chemistry analyzer

CE IVD

The NGAL Test™ Calibrator Kit

REF ST002CA

Revision: TNT2011-10-EN

INTENDED USE

The NGAL Test™ is a particle-enhanced turbidimetric immunoassay for the quantitative determination of neutrophil gelatinase-associated lipocalin (NGAL) in human urine, EDTA plasma and heparin plasma on automated clinical chemistry analyzers. The NGAL Test™ Calibrator Kit contains calibrator solutions for establishing the calibration curve. NGAL measurements are useful in the diagnosis of acute kidney injury which may lead to acute renal failure.

REAGENTS

Composition

Each Calibrator (1-5) contains 1 mL of a ready-to-use solution of recombinant human NGAL in a HEPES buffer containing preservative.

VALUE ASSIGNMENT

No internationally approved reference material for NGAL is currently available. The NGAL solutions in The NGAL Test™ Calibrator Kit have been value-assigned by turbidimetry using a precise transfer protocol ensuring traceability to the BioPorto Diagnostics master calibrator. The reference material was value assigned by measurement of light absorbance at 280 nm using a theoretically calculated extinction coefficient based on the amino acid composition.

WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use IVD

- Read all instructions before starting this test
- This kit should only be used by qualified laboratory staff
- Do not pipette by mouth
- Do not shake the calibrators
- Do not pour calibrators back into their original containers once transferred
- Do not switch caps on containers as it may cause contamination or mix-up
- Do not use after the expiry date on the labels
- All solutions supplied should be handled carefully and disposed of in accordance with national and local regulations
- All equipment used in this test should be sterilized by one of the following methods:
 - Soak in 3.5% v/v glutaraldehyde for 30 minutes or longer
 - Soak in 0.5% w/v sodium hypochlorite for 1 hour or longer
 - Autoclave at 121°C for 20 minutes or longer

STORAGE AND STABILITY

Shelf life at 2-8°C: See expiry date on the label

Stability after opening: 4 weeks at 2-8°C, capped. On the Hitachi 917 analyzer, the manufacturer obtains satisfactory performance for up to 8 weeks.

PROCEDURE

Materials provided

Calibrator 1 (150 ng/mL), 1 mL
Calibrator 2 (600 ng/mL), 1 mL
Calibrator 3 (1500 ng/mL), 1 mL
Calibrator 4 (3000 ng/mL), 1 mL
Calibrator 5 (5000 ng/mL), 1 mL

MATERIALS REQUIRED BUT NOT PROVIDED

- The NGAL Test™ Reagent Kit REF ST001CA
- The NGAL Test™ Control Kit REF ST003CA
- 0.9% w/v aqueous sodium chloride solution as zero calibrator
- Clinical chemistry analyzer
- Analyzer-specific application note (available for a number of analyzers)
- Analyzer-specific reagent containers

CALIBRATION PROCEDURE

Use the Calibrators, with aqueous 0.9% w/v sodium chloride solution as the zero calibrator, to establish a 6-point calibration curve as described in the instrument manual. The assigned NGAL concentration values of the calibrators are given on the labels of the calibrator vials.

Each laboratory should determine the appropriate frequency of calibration. Calibration should be repeated at least once a month or when a new reagent lot is used.

Drip the required volume into a sample cup and analyze in the same way as for specimens.

Measurements should be carried out in accordance with the specific application note for the automated chemistry analyzer used. Please see the Instructions for Use for The NGAL Test™ Reagent Kit REF ST001CA, which also give details of test principle and performance characteristics of The NGAL Test™.

REF Catalogue number

IVD In vitro diagnostic medical device

LOT Batch code

IFU Valid version of instructions

Consult instructions for use

CE European Conformity

Use by

Manufacturer

Temperature limitation

EN

BioPorto Diagnostics A/S
Tuborg Havnevej 15, st.
DK-2900 Hellerup
Denmark

Phone (+45) 4529 0000
Fax (+45) 4529 0001
E-mail info@bioporto.com
Web www.bioporto.com
www.ngal.com

The NGAL Test™

For your clinical chemistry analyzer

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