

**Anti-Glucagon-like peptide-1 (GLP-1, Mid-molecule specific)
 Mouse monoclonal antibody**

Subclass: IgG1/k

PRODUCT NO.

HYB 147-12

Clone: 4F3

PRESENTATION

Preparation: Protein-A/G purified

Content: Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.

Solvent: 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide

Storage: 4-8°C without exposure to light. No precautions necessary during handling.

ANTIGEN

Glucagon-like peptide-1(7-36)amide (GLP-1(7-36)amide) is the principal active form of GLP-1, the other being GLP-1(7-37). GLP-1 is a peptide hormone of the glucagon family, produced by the L cells of the intestinal mucosa from the same prohormone as glucagon. The active forms are potent stimulators of glucose-dependent insulin secretion. The sequence of GLP-1 is fully conserved in all mammalian species examined so far.

IMMUNOGEN

Synthetic GLP-1(7-36)amide coupled to carrier and adsorbed onto aluminum hydroxide gel

SPECIFICITY

Reacts with all forms of GLP-1, including precursor and GLP-1(9-37) /GLP-1(9-36amide) metabolite

EPI TOPE SPECIFICITY

Mid-molecular epitope of GLP-1

REACTIVITY

HYB 147-12 binds to GLP-1 when coated directly onto the microtiter well, and binds GLP-1(7-36)amide in solution giving a K_a of 4.0×10^8 in inhibition ELISA. HYB 147-12 cross-reacts <0.4% with coated glucagon. In inhibition ELISA no binding of free glucagon in solution is detected, giving an estimated cross-reactivity of <0.2%. Although not tested, HYB 147-12 is likely to detect all known molecular forms of GLP-1 in immunohistochemistry.

Biotinylated HYB 147-12 is the preferred detection antibody for measuring C-terminally amidated forms of GLP-1 in combination with HYB 147-06 as capture antibody. HYB 147-12 can be used as a capture antibody in combination with ABS 046-03B as a detection antibody for measuring non-amidated GLP-1 forms and cross-reacting about 16% with C-terminally amidated GLP-1. Results show detection limits of 44pmol/L which is 10-20 times higher than the basal concentration of GLP-1, so the assays have to be optimized.

In Western blotting a dilution guideline of 1/2000 has proved successful.

CULTURE MEDIUM

RPMI 1640 with 10% fetal calf serum

FUSION PARTNER

X63-Ag8.653

IMMUNIZATION

Female BALB/c mice by intraperitoneal injections

APPLICATION

Method	Usability	References
ELISA	Yes	
Immunoblotting	Yes	
Immunohistochemistry	Yes	

REFERENCES

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.