NGAL Rapid ELISA Kit

EN

DE

FR

IT

ES

DA

SE

KIT 037

BIOPORTO[®]
Diagnostics

Please read these instructions carefully**INTENDED USE**

For the *in vitro* determination of human NGAL in urine, plasma or serum as a marker of acute renal injury which may lead to acute renal failure.

INTRODUCTION

NGAL¹ (neutrophil gelatinase-associated lipocalin) belongs to the lipocalin family of proteins. These are typically small secreted proteins characterized by their ability to bind small, hydrophobic molecules in a structurally conserved pocket formed by β -pleated sheet, to bind to specific cell-surface receptors and to form macromolecular complexes. NGAL has many synonyms: it is also known as NL (neutrophil lipocalin); HNL: human NL², lipocalin 2, oncogene protein 24p3³ or uterocalin⁴ (in the mouse) and neuron-related lipocalin⁵ or 25 kDa α_2 -microglobulin-related protein⁶ (in the rat). Human NGAL consists of a single disulfide-bridged polypeptide chain of 178 amino-acid residues with a calculated molecular mass of 22 kDa, but glycosylation increases its apparent molecular mass to 25 kDa. In neutrophils (neutrophilic polymorphonuclear leukocytes) and urine it occurs as monomer, with a small percentage of dimer and trimer, and also in complex with 92-kDa human neutrophil type IV collagenase, also called gelatinase B or matrix metalloproteinase-9 (MMP-9)⁷.

NGAL was originally isolated from the supernatant of activated human neutrophils¹, but it is also expressed at a low level in other human tissues including the kidney, prostate and epithelia of the respiratory and alimentary tracts^{8,9}. It is strongly expressed in adenomas and inflamed epithelia of the bowel¹⁰, adenocarcinomas of the breast¹¹, and urothelial carcinomas¹².

Because of its small molecular size and resistance to degradation, NGAL is readily excreted and detected in the urine, both in its free form and in complex with MMP-9. Urinary levels correlate with plasma or serum levels whatever the cause of increased NGAL production (BioPorto Diagnostics data), but particularly high urinary levels can be

expected when NGAL is released directly into the urine by the kidney tubules or urothelial carcinomas. It is uncertain how far NGAL-MMP-9 complexes from sources remote from the urinary tract are excreted as such into the urine or reform in the urine after independent excretion of NGAL and MMP-9⁷.

While the functions of NGAL are not fully understood, NGAL appears to be upregulated in cells under "stress", e.g. from infection, inflammation, ischemia or neoplastic transformation, or in tissues undergoing involution, such as the postpartum mouse uterus and mammary glands on weaning. In relation to a possible antibacterial role, it binds enterobactin and other siderophores, depriving the microorganisms of Fe³⁺, an important nutritional requirement¹³. Its complex formation with MMP-9 appears to protect MMP-9 enzymatic activity from degradation⁷. The upregulation of NGAL in involuting tissues has led to the postulation of a role in apoptosis, but it appears more likely that NGAL is associated with a survival response¹⁴. This seems to be so in the kidney, where NGAL-siderophore-iron complex rescues the mouse kidney from ischemic injury¹⁵.

NGAL and the kidney. Even before NGAL had been isolated from human neutrophils, its mouse homologue 24p3 was known to be expressed by kidney cells and to undergo an early, dramatic upregulation (14- to 20-fold) in response to SV 40 viral infection¹⁶. A similar early and dramatic upregulation was later observed in rat proximal tubule cells after ischemia-reperfusion injury¹⁷, and raised plasma levels of NGAL were found to be strongly correlated with decreased renal function in patients with renal damage due to systemic vasculitis¹⁸. The results for renal ischemia-reperfusion injury were subsequently confirmed and extended to nephrotoxic agents^{19,20,21}. It has been suggested that urinary NGAL levels may serve as an early marker for ischemic renal injury in children after cardiopulmonary bypass²². Raised urinary and serum NGAL levels have also been observed in patients with established renal failure (BioPorto Diagnostics data) and patients with functioning renal grafts also showed raised urinary levels (detectable by Western

blotting)⁹². It is therefore apparent that a large variety of renal disorders are associated with raised plasma and urinary levels of NGAL. While plasma and urinary NGAL levels are closely correlated in acute conditions, it is to be expected that urinary NGAL levels will be particularly high after ischemic renal injury severe enough to result in acute renal failure, acute tubular necrosis or acute tubulo-interstitial nephropathy. However, the use of urinary NGAL as a marker for these conditions is subject to the proviso that other concurrent conditions that are independently associated with raised NGAL levels must be taken into account.

NGAL in inflammation/infection. NGAL is released from the secondary granules of activated neutrophils¹ and plasma levels rise in inflammatory or infective conditions, especially in bacterial infections⁹³. Thus the level of NGAL in plasma or serum has been proposed as a marker of infection. However, as levels of NGAL may also be raised in neoplastic conditions and renal disorders independently of any infective process, this proposed application should be treated with caution. NGAL may also be raised in infections in patients with an uncountably low number of neutrophils due to leukemia or treated leukemia, showing that the source of the raised NGAL in infections is not only the neutrophils. Indeed, serum NGAL levels correlate very poorly with the neutrophil count in patients with varying degrees of infection or inflammation (BioPorto Diagnostics data).

NGAL and neoplasia. The various types of cancer in which NGAL may be upregulated (often with MMP-9) have been referred to above. This has been shown by its expression in tumor cells and its high urinary levels, both in the free form and complexed with MMP-9⁷. Indeed, it has been proposed that urinary NGAL-MMP-9 complexes may serve as a marker of disease status for breast cancer patients⁹⁴. Plasma levels have not usually been measured in these cases.

PRINCIPLE OF THE ASSAY PROCEDURE

The assay is an ELISA performed in microwells coated with a monoclonal antibody against human NGAL. Bound NGAL is detected with a horseradish peroxidase (HRP)-conjugated monoclonal antibody and the assay is developed by incubation with a color-forming substrate. The assay is a rapid 2-step procedure:

Step 1. Aliquots of calibrators, diluted samples and any controls are incubated with HRP-conjugated detection antibody in the coated microwells. Only NGAL will bind to both coat and detection antibody, while unbound materials are removed by washing.

Step 2. A chromogenic peroxidase substrate containing tetramethylbenzidine (TMB) is added to each test well. The HRP linked to the bound detection antibody reacts with the substrate to generate a colored product. The enzymatic reaction is stopped chemically, and the color intensity is read at 450 nm in an ELISA reader. The color intensity (absorbance) is a function of the concentration of NGAL originally added to each well. The results for the calibrators are used to construct a calibration curve from which the concentrations of NGAL in the test specimens are read.

KIT COMPONENTS

Item	Contents	Quantity
①	12 x 8 coated Microwells + Frame	96 wells
②	5x Sample Diluent Conc.	1 x 60 mL
③-⑤	NGAL Rapid Calibrator 1-6 The exact concentration of each calibrator is printed on the label of the vials and on the QC-certificate.	6 x 1 mL
④	25x Wash Solution Conc.	1 x 30 mL
⑤	HRP-conjugated NGAL Antibody	1 x 6 mL
⑥	TMB Substrate	1 x 12 mL
⑦	Stop Solution	1 x 16 mL
⑧	Polypropylene U-Microwell Plate	96 wells

Note: Liquid reagents contain the preservatives sodium azide, thimerosal or Bronidox L. These may be harmful if ingested.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Adjustable micropipettes covering the range 1-1000 μ L and corresponding disposable pipette tips
2. Polypropylene tubes to contain up to 1000 μ L
3. Tube racks
4. Adjustable 8- or 12-channel micropipette (50-250 μ L range) or repeating micropipette (optional)
5. Clean 1 L and 500 mL graduated cylinders
6. Deionized or distilled water
7. Cover for microplate
8. Clean container for diluted Wash Solution
9. Apparatus for filling wells during washing procedure (optional)
10. Lint-free paper towels or absorbent paper
11. Disposable pipetting reservoirs
12. Timer (60-minute range)
13. Calibrated ELISA plate reader capable of reading at 450 nm (preferably subtracting reference values at 650 or 620 nm)
14. Sodium hypochlorite (household bleach 1:10 dilution) for decontamination of specimens, reagents and materials

PRECAUTIONS

For *in vitro* diagnostic use only

1. This kit should only be used by qualified laboratory staff.
2. Use separate pipette tips for each sample, calibrator and reagent to avoid cross-contamination.
3. Use separate reservoirs for each reagent. This applies especially to the TMB Substrate.
4. After use, decontaminate all specimens, reagents and materials by soaking for at least 30 minutes in sodium hypochlorite solution (household bleach diluted 1:10).
5. To avoid droplet formation during washing, aspirate the wash solution into a bottle containing bleach.
6. Avoid release into the environment. Dispose of containers and unused contents in a safe way and in accordance with national and local regulations.
7. The 5x Sample Diluent Concentrate is

preserved with 0.25% sodium azide (corresponding to 0.05% in final diluted solution). Sodium azide is harmful in contact with skin and if swallowed. Sodium azide is harmful to aquatic organisms and may cause long-term adverse effects in the aquatic environment. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). According to EU regulations, no danger labeling is necessary for the diluted solution.

8. The Stop Solution contains 0.5 mol/L sulfuric acid and can cause irritation or burns to the skin and eyes. If contact occurs, rinse immediately with plenty of water and seek medical advice.
9. Do not interchange components from kits with different batch numbers. The components have been standardized as a unit for a given batch.
10. Hemolyzed, hyperlipemic, heat-treated or contaminated specimens may give erroneous results.
11. Do not dilute clinical specimens directly in the coated microwells.
12. Do not touch or scrape the bottom of the coated microwells when pipetting or aspirating fluid.
13. Incubation times and temperatures other than those specified may give erroneous results.
14. Do not allow the wells to dry once the assay has begun.
15. The TMB Substrate is light sensitive. Keep away from bright light.
16. Do not reuse microwells or pour reagents back into their bottles once dispensed.

STABILITY AND STORAGE

1. Store the kit with all reagents at 2-8°C. Do not freeze.
2. Use all reagents before the expiry date on the vial labels.
3. Diluted Wash Solution Concentrate remains stable for 4 weeks at 2-8°C. If not all wells are to be used, dilute only the portion of Wash Solution Concentrate required.
4. Diluted Sample Diluent Concentrate remains stable for 5 days at 2-8°C. If not all wells are

- to be used, dilute only the portion of Sample Diluent Concentrate required.
- For subsequent use, store unused wells in the foil pouch with the desiccant provided and reseal. Always allow foil pouch to equilibrate to room temperature before opening to avoid condensation in/on the coated microwells.

COLLECTION OF SPECIMENS

Handle and dispose of all blood-derived or urine specimens as if they were potentially infectious.

See Precautions, sections 2, 4 and 5.

Determination of NGAL in a single specimen requires 10 μ L of urine, plasma or serum. Blood specimens should be collected aseptically into EDTA tubes or plain tubes by qualified staff using approved venepuncture techniques. Plasma or serum should be prepared by standard techniques for clinical laboratory testing. Urine should be centrifuged. Cap the prepared clinical specimens. If the assay cannot be performed within 24 hours or specimens are to be shipped, freeze the specimens at -20°C or below. For long-term storage of clinical specimens, -70°C or below is recommended. Avoid repeated freezing and thawing. Do not use hemolyzed, hyperlipemic, heat-treated or contaminated specimens.

PREPARATION OF REAGENTS

- Bring all specimens and reagents to room temperature (20 - 25°C). Mix specimens thoroughly by gentle inversion and if necessary clear visible particulate matter by low-speed centrifugation.
- Determine the number of specimens to be tested (in duplicate) plus any internal laboratory control specimens (in duplicate) plus any reagent blank wells. The precoated wells can be used as strips of 8 or as individual wells. Single wells are handled by breaking individual wells apart and placing each well in the frame at an appropriate position. Letters and notches on the wells allow the individual wells to be identified. Add 12 wells for the 6 calibrators (in duplicate). Remove the number of microwells required and replace the remainder in the foil pouch with desiccant at 2 - 8°C .


- Wash Solution: Dilute the 25x Wash Solution Concentrate by pouring the total contents of the bottle (30 mL) into a 1-L graduated cylinder and add distilled or deionized water to a final volume of 750 mL. Mix thoroughly and store at 2 - 8°C after use. If not all the wells are to be used, dilute only the required volume of Wash Solution Concentrate with 24 volumes of water to produce a 1/25 dilution.
- Sample Diluent: Dilute the 5x Sample Diluent Concentrate (contains yellow dye to aid pipetting) by pouring the total contents of the bottle (60 mL) into a 500-mL graduated cylinder and add distilled or deionized water to a final volume of 300 mL. Mix thoroughly and store at 2 - 8°C after use. If not all the wells are to be used, dilute only the required volume of Sample Diluent Concentrate with 4 volumes of water to produce a 1/5 dilution.
- NGAL Rapid Calibrators (ready to use): The assigned concentrations are indicated on their labels. Do not dilute further.
- HRP-conjugated NGAL Antibody (ready to use): Do not dilute further.
- TMB Substrate (ready to use): Do not dilute further.
- Stop Solution (ready to use): Do not dilute further.

PREPARATION OF SAMPLES

The approximate range of the standard curve is 0.2-20 ng/mL and the diagnostically relevant range is 100-500 ng/mL for plasma or serum and 50-500 ng/mL for urine (see Interpretation of Results). Therefore an initial screening dilution of 1/100 for plasma or serum and 1/50 for urine is recommended. The 1/100 dilution can be prepared by diluting 10 μ L of plasma or serum in 990 μ L of Sample Diluent and the 1/50 dilution can be prepared by diluting 10 μ L of urine in 490 μ L of Sample Diluent. Dilutions are mixed by inversion or moderate vortexing. Reassay of out-of-range samples at lower or higher dilution may be necessary. Dilutions lower than 1/25 should not be used.

ASSAY PROCEDURE**(See also schematic overview)**

1. Prepare the assay protocol, assigning the appropriate wells for setting up calibrators, diluted patient specimens and any internal laboratory controls in duplicate. If a reference wavelength of 650 or 620 nm is not available on the ELISA reader, a reagent blank well can be assigned. This is set up with 50 μ L of Sample Diluent instead of diluted sample and processed like the other wells.
2. Dilute samples according to the expected NGAL concentrations (1/100 for plasma or serum and 1/50 for urine will be suitable for most samples).
3. Pipette a sufficient volume of each calibrator, each diluted sample and any internal laboratory controls into the appropriate wells of the polypropylene U-microwell plate to permit subsequent transfer of 50 μ L volumes to corresponding coated microwells.
4. Pipette 50 μ L volumes of HRP-conjugated NGAL Antibody into the corresponding positions in the coated microwells. Then with a multichannel pipette rapidly transfer 50 μ L volumes of the calibrator solutions, diluted samples and internal controls from the U-wells into the corresponding coated wells already containing the detection antibody. This method of sample addition is recommended to reduce the difference in incubation time between the first and last samples added to the coated microwells.
5. Cover the wells and incubate for **30 minutes** at room temperature on a shaking platform set at 200/minute.
6. Aspirate the contents of the microwells and wash the microwells three times with at least 300 μ L of the previously diluted Wash Solution. If washing is performed manually, empty the microwells by inversion and gentle shaking into a suitable container, followed by blotting in the inverted position on a paper towel. A dwell time of 1 minute before emptying is recommended for at least the last wash of the cycle. The vigor with which diluted Wash Solution is filled into or emptied from the wells influences final color development. Manual pipetting, which may be very gentle and lead to high color development, is only recommended in the absence of alternatives such as filling the wells by immersion, using a multi-channel manual washing dispenser, or using an automatic washing apparatus.
7. Dispense 100 μ L of TMB Substrate (ready to use) into each microwell. The use of a multichannel micropipette is recommended to reduce pipetting time. Cover the wells and incubate for **exactly 15 minutes** at room temperature in the dark. Start the clock when filling the first well.
8. Add 100 μ L of Stop Solution (ready to use) to each well, maintaining the same pipetting sequence and rate as in Step 7. Mix by gentle shaking for 20 seconds, avoiding splashing. Read the wells within 30 minutes.
9. Read the absorbances of the wells at 450 nm in an appropriate microplate reader (reference wave-length 650 or 620 nm). If no reference wavelength is available, the value of the reagent blank well is subtracted from each of the other values before other calculations are performed.

Schematic overview of assay procedure**CALCULATION OF RESULTS**

For calculation of results the exact values of the calibrators printed on the labels of the calibrators and on the QC-certificate should be used.

The basic principle is to construct a calibration curve by plotting the mean of duplicate absorbance values for each NGAL Rapid Calibrator on the y-axis against the corresponding NGAL concentrations in ng/mL on the x-axis. The calibration curve must meet the validation requirements. The NGAL concentration of each diluted sample is then found by locating the point on the curve corresponding to the mean of

duplicate absorbance values for the diluted sample and reading its corresponding concentration in ng/mL from the x-axis. The concentration of NGAL in the undiluted specimen is calculated by multiplying this result by the sample dilution factor.

This procedure can be performed manually using graph paper with linear x and y scales. A smooth curve can be drawn through the points or adjacent points can be joined by straight lines. The latter procedure may slightly overestimate concentration values between points when the curve is slightly convex to the left, which is the typical finding. Although the curve may approximate to a straight line, it is both practically and theoretically incorrect to calculate and draw the straight line of best fit and to read the results from this.

The procedure can also be performed by an ELISA reader software program incorporating curve fitting procedures. The procedure of choice is to use linear x and y axes with 4-parameter logistic curve fitting. Diluted samples that give a mean absorbance above that for the NGAL Rapid Calibrator 6 or below that for the NGAL Rapid Calibrator 2 are out of the range of the assay and their concentrations should be noted as >A ng/mL and <B ng/mL, respectively, where A is the exact concentration of NGAL Rapid Calibrator 6 and B is the exact concentration of NGAL Rapid Calibrator 2. The corresponding concentrations in the undiluted sera are calculated >(A x dilution factor) ng/mL and <(B x dilution factor) ng/mL, respectively. If necessary, these samples can be re-assayed at higher and lower dilutions for high- and low-reading samples, respectively. The new dilution factors should be those estimated to give absorbance values that fall well within the range of the calibration curve, but dilutions lower than 1/25 should not be used.

VALIDATION OF CALIBRATION CURVE

The mean absorbance for the NGAL Rapid Calibrator 6 should be >1.5. The mean absorbance for any NGAL Rapid Calibrator should be higher than that for the previous NGAL Rapid Calibrator, e.g. absorbance(NGAL Rapid Calibrator 6) > absorbance(NGAL Rapid Calibrator 5). The curve

should be slightly convex to the left when the results are plotted on linear axes.

Out-of-line points for individual calibrators: One or more individual calibrators may give anomalous absorbance readings. One or both of the duplicate values may be out of line, and the mean of the duplicates may be out of line. This error is significant if it impairs satisfactory curve fitting by the 4-parameter logistic method, which, as a result of the anomalous value, is shifted away from other calibrator points that are in fact correct. The calibrator points and fitted curve should always be examined for correct fit before any calculations of concentration from it are accepted. A poorly fitting curve will also be revealed by a high sum of residual squares. If only one calibrator is affected, which is not the highest calibrator, two courses of action are possible:

i) An erroneous singlet or duplicate result should be eliminated from the curve, and the remaining results refitted by the 4-parameter logistic procedure. If a satisfactory fit is obtained, provisional concentration results can be calculated from it.

ii) If no satisfactory fit can be obtained in this way, but the curve is otherwise consistent, provisional results can be obtained from straight lines or simple cubic spline fitting between the means of duplicates, omitting the erroneous point.

If two or more calibrators are affected, the assay should be repeated.

A deviant result for an individual calibrator can be due to operator error or to calibrator deterioration. If both duplicate values are consistently out of line in successive assays, the calibrator is faulty and should be omitted.

TRACEABILITY OF CALIBRATOR VALUE

The NGAL concentration of diluted calibrator material has been assigned by comparison in ELISA with a purified preparation of recombinant human NGAL quantified by Rigshospitalet, Copenhagen, Denmark.

INTERPRETATION OF RESULTS

The finding of a raised urinary, plasma or serum

level of NGAL cannot be independently diagnostic of any single pathology. As stated in the Introduction, a variety of independent pathologies are associated with raised levels of urinary or plasma NGAL. Physicians must interpret the significance of any raised NGAL level in the light of each patient's clinical features.

The capacity of urinary and plasma NGAL levels to diagnose acute renal failure presumed to result from renal injury was determined in a series of 60 consecutive patients admitted to a non-specialized intensive care unit and classifiable by other criteria as showing or not showing signs of renal failure. Many of these patients had concurrent sepsis and some had cancers. Raised maximal NGAL levels were associated with acute renal failure as shown:

Acute renal failure

Sample	Urine	Plasma
Cutoff value	350 ng/mL	400 ng/mL
Diagnostic specificity	89.3%	96.3%
Diagnostic sensitivity	96.9%	84.8%
Positive predictive value	91.2%	93.1%
Negative predictive value	96.2%	83.9%

An NGAL concentration in urine above 350 ng/mL or in plasma above 400 ng/mL is highly likely (>90%, positive predictive value) to be due to renal injury that may lead to acute renal failure. NGAL concentrations below the cutoff values may indicate either no or only slight renal injury in this patient population. A sudden rise in the urinary NGAL concentration above the patient's previous values may indicate renal injury even if the cutoff value is not exceeded, and close surveillance of renal function is recommended.

QUALITY CONTROL

Laboratories intending to perform repeated assays should establish their own high-reading and low-reading control sera or urine, stored in small (e.g. 50 µL) aliquots at -70°C or below. An aliquot of each should be thawed and tested in each assay and a record kept of successive results. This serves as a control of test performance, test integrity and operator reliability. The results should be examined for drift (tendency for successive results to rise or fall) or significant deviation from the mean of previous results. Values not deviating by more than 20% from the mean of previous values can be taken to indicate acceptability of the assay. Aliquots of control serum should not be refrozen for repeated assay once thawed, and if a further assay is performed, fresh control aliquots and fresh dilutions of patient specimens should be used.

LIMITATIONS

The finding of a raised urinary, plasma or serum level of NGAL cannot be independently diagnostic of any single pathology. As stated in the Introduction, a variety of independent pathologies are associated with raised levels of urinary or plasma NGAL. Physicians must interpret the significance of any raised NGAL level in the light of each patient's clinical features.

The cutoff values given in the Interpretation of Results are based on the maximal values in the patient population described, which was characterized by the presence of multiple confounding factors. Other cutoffs, perhaps lower, may be appropriate for other patient groups with fewer confounding factors. Sudden rises in NGAL may indicate acute renal injury even if they do not exceed the cutoff values.

EXPECTED RESULTS

The mean NGAL concentration in samples from healthy donors was 63 ng/mL (range 37–106 ng/mL, n = 80) in EDTA plasma and 5.3 ng/mL (range 0.7–9.8 ng/mL, n = 7) in urine. In unselected patients admitted to intensive care, the NGAL concentrations in urine ranged from 9 ng/mL to 40,000 ng/mL

(40 µg/mL) in urine (n = 60) and from 25 ng/mL to 3491 ng/mL in EDTA plasma (n = 60).

PERFORMANCE CHARACTERISTICS

Limit of detection: The lowest concentration of NGAL giving an absorbance reading greater than 2 SD above the mean zero (NGAL Rapid Calibrator 1) reading (n = 6) was determined 4 times by different operators. The results were <0.1 ng/mL, these being lower than the value of NGAL Rapid Calibrator 2.

Precision: Intraassay variation was determined by measurement of NGAL in a urine sample and a serum sample with 6 replicates in 4 separate assays performed by 4 different operators. The following results were obtained (CV = coefficient of variation):

Samples	CV Median (range)
Urine	3.4% (2.5-8.4)
Serum	3.2% (1.5-5.4)

Interassay variation was determined by measurement of NGAL in 7 diluted urine samples, 6 diluted serum and 7 diluted EDTA plasma samples with 2 replicates in 4 separate assays, performed by 4 different operators. The following results were obtained:

Samples	CV Median (range)
Urine	12.6% (5.7-30.4)
Serum	7.4% (4.9-10.5)
EDTA plasma	13.5% (6.5-22.3)

Analytical recovery: Urine and plasma samples were spiked with recombinant human NGAL and analyzed in the assay. Recovery was calculated from (Measured/Expected) expressed as a percentage.

Sample	Measured	Expected	Recovery
Urine	3.4 ng/mL	3.5 ng/mL	96%
Plasma	5.4 ng/mL	6.2 ng/mL	85%

Linearity: NGAL was measured in serial dilutions ($n = 6$) of a urine sample and a serum sample in 4 separate assays performed by 4 different operators. The CV of the mean of the measured values corrected for the dilution was 7.0% for urine and 8.5% for serum, demonstrating satisfactory linearity.

Sample material: Analysis of samples of serum, EDTA plasma or urine showed no significant differences in analytical recovery, linearity or precision. However, NGAL concentrations are often slightly higher in serum than in plasma²³, probably because of the release of NGAL from neutrophils during blood clotting.

Specificity: The two monoclonal antibodies against human NGAL used in the assay have been shown to bind to different preparations of recombinant human NGAL and to give a single band at 25 kDa on Western blot analysis of a reduced postnuclear supernatant from human neutrophils²⁵.

LIABILITY

For *in vitro* diagnostic use in the European Union and Canada only. For research use only in the rest of the world. Approval for diagnostic use in the rest of the world has not been sought.

This kit is only intended for the *in vitro* determination of NGAL in human plasma, serum or urine.

The kit is only intended for use by qualified personnel carrying out research or diagnostic activities.

If the recipient of this kit passes it on in any way to a third party, this instruction must be enclosed, and said recipient shall at own risk secure in favor of BioPorto Diagnostics A/S all limitations of liability herein.

BioPorto Diagnostics A/S shall not be responsible for any damages or losses due to using the kit in any way other than as expressly stated in these Instructions.

The liability of BioPorto Diagnostics A/S shall in no event exceed the commercial value of the kit.

BioPorto Diagnostics A/S shall under no

circumstances be liable for indirect, special or consequential damages, including but not limited to loss of profit.

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1. Kjeldsen L, Johnsen AH, Sengelov H, Borregaard N (1993) Isolation and primary structure of NGAL, a novel protein associated with human neutrophil gelatinase. *J Biol Chem* 268:10425-10432.
2. Xu SY, Carlson M, Engstrom A, Garcia R, Peterson CG, Venge P (1994) Purification and characterization of a human neutrophil lipocalin (HNL) from the secondary granules of human neutrophils. *Scand J Clin Lab Invest* 54:365-376.
3. Flower DR, North AC, Attwood TK (1991) Mouse oncogene protein 24p3 is a member of the lipocalin protein family. *Biochem Biophys Res Commun* 180:69-74.
4. Liu Q, Ryon J, Nilsen-Hamilton M (1997) Uterocalin: a mouse acute phase protein expressed in the uterus around birth. *Mol Reprod Dev* 46:507-514.
5. Stoesz SP, Gould MN (1995) Overexpression of neu-related lipocalin (NRL) in neu-initiated but not ras or chemically initiated rat mammary carcinomas. *Oncogene* 11:2233-2241.
6. Triebel S, Blaser J, Reinke H, Tschesche H (1992) A 25 kDa alpha 2-microglobulin-related protein is a component of the 125 kDa form of human gelatinase. *FEBS Lett* 314:386-388.
7. Yan L, Borregaard N, Kjeldsen L, Moses MA (2001) The high molecular weight urinary matrix metalloproteinase (MMP) activity is a complex of gelatinase B/MMP-9 and neutrophil gelatinase-associated lipocalin (NGAL). Modulation of MMP-9 activity by NGAL. *J Biol Chem* 276:37258-37265.
8. Cowland JB, Borregaard N (1997) Molecular characterization and pattern of tissue expression of the gene for neutrophil gelatinase-associated lipocalin from humans. *Genomics* 45:17-23.
9. Friedl A, Stoesz SP, Buckley P, Gould MN (1999) Neutrophil gelatinase-associated lipocalin in normal and neoplastic human tissues. Cell type-specific pattern of expression. *Histochem J* 31:433-441.
10. Nielsen BS, Borregaard N, Bundgaard JR, Timshel S, Sehested M, Kjeldsen L (1996) Induction of NGAL synthesis in epithelial cells of human colorectal neoplasia and inflammatory bowel diseases. *Gut* 38:414-420.
11. Stoesz SP, Friedl A, Haag JD, Lindstrom MJ, Clark GM, Gould MN (1998) Heterogeneous expression of the lipocalin NGAL in primary breast cancers. *Int J Cancer* 79:565-572.
12. Monier F, Surla A, Guillot M, Morel F (2000) Gelatinase isoforms in urine from bladder cancer patients. *Clin Chim Acta* 299:11-23.
13. Goetz DH, Holmes MA, Borregaard N, Bluhm ME, Raymond KN, Strong RK (2002) The neutrophil lipocalin NGAL is a bacteriostatic agent that interferes with siderophore-mediated iron acquisition. *Mol Cell* 10:1033-1043.
14. Tong Z, Wu X, Ovcharenko D, Zhu J, Chen CS, Kehrer JP (2005) Neutrophil gelatinase associated lipocalin as a survival factor. *Biochem J* 391:441-448.
15. Mori K, Lee HT, Rapoport D, Drexler IR, Foster K, Yang J, Schmidt-Ott KM, Chen X, Li JY, Weiss S, Mishra J, Cheema FH, Markowitz G, Suganami T, Sawai K, Mukoyama M, Kunis C, D'Agati V, Devarajan P, Barasch J (2005) Endocytic delivery of lipocalin-siderophore-iron complex rescues the kidney from ischemia-reperfusion injury. *J Clin Invest* 115:610-621.
16. Hrabá-Renevey S, Turler H, Kress M, Salomon C, Weil R (1989) SV40-induced expression of mouse gene 24p3 involves a post-transcriptional mechanism. *Oncogene* 4:601-608.
17. Matthaeus T, Schulze-Lohoff E, Ichimura T, Weber M, Andreucci M, Park KM, Alessandrini A, Bonventre JV (2001) Co-regulation of neutrophil gelatinase-associated lipocalin and matrix metalloproteinase-9 in the posts ischemic

- rat kidney. *J Am Soc Nephrol* 12:787A.
18. Ohlsson S, Wieslander J, Segelmark M (2003) Increased circulating levels of proteinase 3 in patients with anti-neutrophilic cytoplasmic auto-antibodies-associated systemic vasculitis in remission. *Clin Exp Immunol* 131:528-535.
 19. Mishra J, Ma Q, Prada A, Mitsnefes M, Zahedi K, Yang J, Barasch J, Devarajan P (2003) Identification of neutrophil gelatinase-associated lipocalin as a novel early urinary biomarker for ischemic renal injury. *J Am Soc Nephrol* 14:2534-2543.
 20. Amin RP, Vickers AE, Sistare F, Thompson KL, Roman RJ, Lawton M, Kramer J, Hamadeh HK, Collins J, Grissom S, Bennett L, Tucker CJ, Wild S, Kind C, Oreffo V, Davis JW 2nd, Curtiss S, Naciff JM, Cunningham M, Tennant R, Stevens J, Car B, Bertram TA, Afshari CA (2004) Identification of putative gene based markers of renal toxicity. *Environ Health Perspect* 112:465-479.
 21. Mishra J, Mori K, Ma Q, Kelly C, Barasch J, Devarajan P (2004) Neutrophil gelatinase-associated lipocalin: a novel early urinary biomarker for cisplatin nephrotoxicity. *Am J Nephrol* 24:307-315.
 22. Mishra J, Dent C, Tarabishi R, Mitsnefes MM, Ma Q, Kelly C, Ruff SM, Zahedi K, Shao M, Bean J, Mori K, Barasch J, Devarajan P (2005) Neutrophil gelatinase-associated lipocalin (NGAL) as a bio-marker for acute renal injury after cardiac surgery. *Lancet* 365(9466):1231-1238.
 23. Xu SY, Pauksen K, Venge P (1995) Serum measurements of human neutrophil lipocalin (HNL) discriminate between acute bacterial and viral infections. *Scand J Clin Lab Invest* 55:125-131.
 24. Fernandez CA, Yan L, Louis G, Yang J, Kutok JL, Moses MA (2005) The matrix metalloproteinase-9/neutrophil gelatinase-associated lipocalin complex plays a role in breast tumor growth and is present in the urine of breast cancer patients. *Clin Cancer Res* 11:5390-5395.
 25. Kjeldsen L, Koch C, Arnljots K, Borregaard N (1996) Characterization of two ELISAs for NGAL, a newly described lipocalin in human neutrophils. *J Immunol Methods* 196:155-164.

REF

Catalogue number
Bestellnummer
Réf rence du catalogue
Numero di catalogo
N mero de cat logo
Katalognummer
Katalognummer

Use by
Verwendbar bis
Utiliser jusque
Utilizzare entro
Fecha de caducidad
Anvend for
Anv nd f re

IVD

In vitro diagnostic medical device
In-Vitro-Diagnostikum
Dispositif m dical de diagnostic in vitro
Dispositivo medico-diagnostico in vitro
Producto sanitario para diagn stico in vitro
Medicinsk udstyr til in vitro-diagnostik
Medicintekniske produkter f r in vitro diagnostik

Manufacturer
Hersteller
Fabricant
Fabbricante
Fabricante
Fabrikant
Tillverkare

LOT

Batch code
Chargenbezeichnung
Code du lot
Codice del lotto
C digo de lote
Lotnummer
Lot nummer

Keep away from sunlight
Vor Sonnenlicht sch tzen
Ne pas exposer aux rayons solaires
Evitare l'esposizione alla luce del sole
Mantener fuera de la luz solar
M  ikke uds ttes for sollys
Uts tts ej f r direkt solljus

Consult instructions for use
Gebrauchsanweisung beachten
Consulter les instructions d'utilisation
Consultare le istruzioni per l'uso
Consulte las instrucciones de uso
Se brugsanvisning
Se handhavandebeskrivningen

Temperature limitation
Temperaturbegrenzung
Limites de temp rature
Limiti di temperatura
L mite de temperatura
Temperaturbegr nsning
Temperaturbegr nsning

European Conformity
CE-Konformit tskennzeichnung
Conformit  aux normes europ ennes
Conformit  europea
Conformidad europea
Europ isk overensstemmelse
Europeisk  versensst mmelse

Do not reuse
Nicht zur wiederverwendung
Ne pas r utiliser
Non riutilizzare
No reutilizar
M  ikke genbruges
 teranv nd ej

Caution, consult accompanying documents
 Achtung, Begleitdokumente beachten
 Attention, voir notice d'instructions
 Attenzione, vedere le istruzioni per l'uso
 Atención, ver instrucciones de uso
 Forsigtig, se brugsanvisning
 Försiktighet, se handhavandebeskrivningen

Biological risk
 Biogefährdung
 Risques biologiques
 Rischio biologico
 Riesgo biológico
 Biologisk fare
 Biologisk risk

Do not use if package is damaged
 Inhalt beschädigter Packung nicht verwenden
 Ne pas utiliser si l'emballage est endommagé
 Non utilizzare se la confezione è danneggiata
 No usar si el paquete está dañado
 Må ikke anvendes hvis emballagen er beskadiget
 Använd inte om förpackningen är skadad

WASH SOLUTION 25X

Concentrated Wash Solution. Dilute before use.
 Konzentrierte Waschlösung. Vor Gebrauch verdünnen.
 Solution de bain concentrée. Diluer avant usage.
 Soluzione di lavaggio concentrata Diluire prima dell'uso.
 Tampón de lavado concentrado. Diluir antes de usar.
 Vaskeopløsningskoncentrat. Fortyndes før anvendelse.
 Koncentrerad tvättlösning. Späda före användning.

SAMPLE DILUENT 5X

Concentrated Sample Diluent. Dilute before use.
 Konzentriertes Probenverdünnungsmittel. Vor Gebrauch verdünnen.
 Diluant échantillon concentrée. Diluer avant usage.
 Diluente per campione concentrata. Diluire prima dell'uso.
 Tampón de dilución concentrado. Diluir antes de usar.
 Prøvefortyndingskoncentrat. Fortyndes før anvendelse.
 Koncentrad provspädningsvätska. Späda före användning.

BIOPORTO[®]
Diagnostics

BioPorto Diagnostics A/S
Grusbakken 8
DK-2820 Gentofte
Denmark

Phone (+45) 4529 0000
Fax (+45) 4529 0001
E-mail info@bioporto.com
Web www.bioporto.com