

Anti-L-Ficolin (human)**Mouse monoclonal antibody**

Subclass: IgG2a/k

CAT. NO.

ABS 005-19

Clone:19

SPECIFICITY	ABS 005-19 reacts specifically with L-Ficolin. No cross reactivity is seen with M-Ficolin (Ficolin-1) and H-Ficolin (Ficolin-3).
IMMUNOGEN	Purified procaryot expressed recombinant human L-Ficolin and full-length L-Ficolin (Ficolin-2) protein
TESTED APPLICATIONS	ELISA, WB
SPECIES REACTIVITY (POSITIVE)	Human
SPECIES REACTIVITY (NEGATIVE)	Not determined
EPITOPE SPECIFICITY	Not determined

PRESENTATION

Content:	Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.
Preparation:	Protein-A purified
Form:	Liquid
Solvent:	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
Storage:	4-8°C without exposure to light. No precautions necessary during handling.

APPLICATION	ELISA: ABS 005-19 is especially well suited as detecting antibody in a sandwich ELISA in combination with ABS 005-16 (1). Furthermore, in ELISA, ABS 005-19 reacts with recombinant L-Ficolin coated directly into the microtiter well. WB: ABS 005-19 can be used in Western blotting (1, 2, 3). IP: ABS 005-19 can be used in immunoprecipitation (4).
-------------	---

TARGET	L-Ficolin (Ficolin-2) is an innate immunity pattern recognition molecule found in human serum. L-Ficolin binds to distinct pathogen-associated molecular patterns (PAMP), such as carbohydrates (GlcNAc), lipoteichoic acid and acetylated groups thereby facilitating phagocytosis and activation of complement through the lectin route using the same serine proteases as mannose-binding lectin (MBL) named MASP ((MBL)-associated serine protease) (1).
--------	--

REFERENCES	<ol style="list-style-type: none">1. Munthe-Fog L, Hummelshøj T, Hansen BE, Koch C, Madsen HO, Skjødt K, Garred P (2007) The impact of FCN2 polymorphisms and haplotypes on the Ficolin-2 serum levels. Scand J Immunol 65:383-392.2. Jensen ML, Honoré C, Hummelshøj T, Hansen BE, Madsen HO, Garred P (2006) Ficolin-2 recognizes DNA and participates in the clearance of dying host cells. Mol Immunol 44:856-865.3. Hummelshøj T, Thielens NM, Madsen HO, Arlaud GJ, Sim RB, Garred P (2006) Molecular organization of human Ficolin-2. Mol Immunol 44:401-411.4. Skjøedt MO, Hummelshøj T, Palarasah Y, Honoré C, Koch C, Skjødt K, Garred P (2010) A Novel Mannose-binding Lectin/Ficolin-associated Protein Is Highly Expressed in Heart and Skeletal Muscle Tissues and Inhibits Complement Activation. J Biol Chem 285:8234-8243.
------------	--

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.