

Anti-D-dimer (human, dog)
Mouse monoclonal antibody
ABS 015-22

Subclass: IgG1/k
Clone: 22

CAT. NO.

SPECIFICITY	ABS 015-22 binds human D-dimer. Cross-reacts < 2% with directly coated D-monomer but not with directly coated fibrinogen.
IMMUNOGEN	Purified human D-dimer
TESTED APPLICATIONS	ELISA, WB
SPECIES REACTIVITY (POSITIVE)	Human, dog
SPECIES REACTIVITY (NEGATIVE)	Not determined
EPITOPE SPECIFICITY	When tested in a pair with itself, with ABS 015-28 and with a fibrinogen-specific polyclonal antibody, the epitope appears to relate to the dimenzation site of D-dimer.

PRESENTATION

Content:	Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.
Preparation:	Protein-A purified
Form:	Liquid
Solvent:	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
Storage:	4-8°C without exposure to light. No precautions necessary during handling.

APPLICATION

ELISA: ABS 015-22 binds free D-dimer in solution.
ABS 015-022 works as capture antibody in sandwich ELISA with biotinylated ABS 015-22 or ABS 015-28 as detection antibody.
Biotinylated ABS 015-22 works as detection antibody in sandwich ELISA for D-dimer with itself (ABS 015-22) or ABS 015-28 as capture antibody.
The combination with ABS 015-28 as capture antibody and ABS 015-22B as detection antibody is recommended.
WB: Due to reducing conditions in Western blotting ABS 015-22 binds the alpha-chain (12kDa) of D-Dimer, and not the entire D-Dimer complex.

TARGET

D-dimer is a fibrin degradation product, a small protein fragment present in the blood after a blood clot is degraded by fibrinolysis.

REFERENCES

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.