

Anti-Procollagen type I C-terminal propeptide (human, PICP)**Mouse monoclonal antibody**Subclass: IgG₁/κ

CAT. NO.

BTE 002-04

Clone: 4

SPECIFICITY	BTE 002-04 is specific for native PICP from amniotic fluid.
IMMUNOGEN	PICP purified from second trimester amniotic fluid by affinity chromatography
TESTED APPLICATIONS	ELISA
SPECIES REACTIVITY (POSITIVE)	Human
SPECIES REACTIVITY (NEGATIVE)	Not determined
EPITOPE SPECIFICITY	Not determined

PRESENTATION

Content:	Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.
Preparation:	Protein-A purified
Form:	Liquid
Solvent:	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
Storage:	4-8°C without exposure to light. No precautions necessary during handling.

APPLICATION **ELISA:** BTE 002-04 reacts with PICP in ELISA.**TARGET** Collagen type I is present in soft connective tissue and bone where it represents approximately 90% of the organic matrix. Collagen type I is produced from procollagen when two teloterminal propeptides are digested from the precursor and released into the circulation. The propeptide levels can be measured in serum as indicators of collagen production. Elevated values of these propeptides in serum are found in patients with various bone metabolism defects. Procollagen type I C-terminal propeptide (PICP) has a molecular weight of approximately 100 kDa and the normal range in serum is 100-300µg/ml.**REFERENCES****CONDITIONS**

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.