

Anti-Tetranectin (human)
Mouse monoclonal antibody
HYB 130-11

Subclass: IgG2a/k
Clone: 5B7

CAT. NO.

SPECIFICITY	HYB 130-11 is specific for amino acids 17-181 of human tetranectin monomer.
IMMUNOGEN	Tetranectin purified from human citrate plasma. Boosted before fusion with recombinant tetranectin produced in <i>E. coli</i> .
TESTED APPLICATIONS	WB, IHC-P, IHC-F, M
SPECIES REACTIVITY (POSITIVE)	Human, bovine
SPECIES REACTIVITY (NEGATIVE)	Not determined
EPITOPE SPECIFICITY	

PRESENTATION

Content:	Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.
Preparation:	Protein-A purified
Form:	Liquid
Solvent:	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
Storage:	4-8°C without exposure to light. No precautions necessary during handling.

APPLICATION

ELISA: HYB 130-11 was used in ELISA
WB: In Western blotting HYB 130-11 reacts strongly with tetranectin monomer (1).
IHC: Strong recognition of tetranectin in immunohistochemistry on both fresh frozen and paraffin-embedded tissues (1, 2).
Microscopy: HYB 130-11 can be used in confocal microscopy and light microscopic immunofluorescence (2).

TARGET

Tetranectin (TN) is a serum and tissue protein, a C-type lectin, which binds to Ca⁺⁺. It is a homotrimer of monomers each with a mass of 20 kDa, plasma or serum concentrations of TN are found to be approximately 10 mg/l (1). In vitro, TN can bind to kringle 4 of plasminogen and enhance the activation of plasminogen to plasmin, catalyzed by tissue plasminogen activator in the presence of poly-D-lysine. TN is best known as a prognostic marker in ovarian cancer.

REFERENCES

- Hogdall CK, Christiansen M, Christensen L, Yazova AK, Koch C, Clemmensen I, Norgaard-Pedersen B (1997) Monoclonal antibodies against human tetranectin, epitope characterization and use in immunohistochemistry. *Clin Chim Acta* 258:159-177.
- Wernbom M, Paulsen G, Nilsen TS, Hisdal J, Raastad T (2012) Contractile function and sarcolemmal permeability after acute low-load resistance exercise with blood flow restriction. *Eur J Appl Physiol* 112:2051-2063.

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.