

## 30 mL 25x Wash Solution Concentrate for KIT 029

PRODUCT NO.	<b>WASH 029-30</b>	
PRESENTATION	Preparation:	25x Wash Solution Concentrate
	Content:	30 mL
	Storage:	Store at 2-8°C. Do not freeze.
	Stability, undiluted:	See label
	Stability after dilution:	4 weeks at 2-8°C. Do not freeze.
PREPARATION	Dilute the 25x Wash Solution Concentrate by pouring the total contents of the bottle (30 mL) into a 1-L graduated cylinder and add distilled or deionized water to a final volume of 750 mL. Mix thoroughly and store at 2-8°C after use.	
INTENDED USE	The 30 mL 25x Wash Solution Concentrate for KIT 029 is intended for use with the MBL Oligomer ELISA Kit when used with automated ELISA washing systems. The kit already contains 30 mL 25x Wash Solution Concentrate, which should be sufficient for manual washing.	
INSTRUCTIONS	Please refer to the instructions for use for KIT 029 MBL Oligomer ELISA Kit for further information.	

### CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.