

**Anti-Glucagon-like peptide-1 (GLP-1(7-36)amide, C-terminal specific)
Mouse monoclonal antibody**

Subclass: IgG1/k

PRODUCT NO.

HYB 147-06

Clone: 8G9

PRESENTATION

Preparation: Protein-A/G purified

Content: Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.

Solvent: 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide

Storage: 4-8°C without exposure to light. No precautions necessary during handling.

ANTIGEN

Glucagon-like peptide 1(7-36)amide (GLP-1(7-36)amide) is the principal active form of GLP-1, the other being GLP-1(7-37). GLP-1 is a peptide hormone of the glucagon family, produced by the L cells of the intestinal mucosa from the same prohormone as glucagon. The active forms are potent stimulators of glucose-dependent insulin secretion. The sequence of GLP-1 is fully conserved in all mammalian species examined so far.

IMMUNOGEN

Synthetic GLP-1(7-36)amide coupled to carrier and adsorbed onto aluminum hydroxide gel

SPECIFICITY

HYB 147-06 is specific for the amidated C-terminus of the peptide and does not react with GLP-1(7-37) (1).

EPI TOPE SPECIFICITY

C-terminal epitope of GLP-1(7-36)amide

REACTIVITY

HYB 147-06 reacts with the amidated C-terminus of GLP-1(7-36)amide, GLP-1(9-36)amide and GLP-1(1-36)amide. HYB 147-06 can be used as capture antibody in sandwich ELISA (1) using HYB 147-12B or ABS 033-10B as detection antibody. HYB 147-06 has been used for the immunoblockade of endogenous GLP-1 in rats (2).

CULTURE MEDIUM

RPMI 1640 with 10% fetal calf serum

FUSION PARTNER

X63-Ag8.653

IMMUNIZATION

Female BALB/c mice immunized by intraperitoneal injections

APPLICATION

Method	Usability	References
ELISA	Yes	1
Immunoblotting	Not determined	
Immunohistochemistry	Yes	

REFERENCES

1. Ghiglione M, Uttenthal LO, Koch C (1993) Monoclonal antibodies to glucagon-like peptide-1. Digestion 54:396-397.
2. van Delft J, Uttenthal O, Koch C, Ghiglione M (1999) Immunoblockade of endogenous glucagon-like peptide-1 by monoclonal antibodies in conscious rats: effect on the insulin response to intragastric glucose. Metabolism 48:41-46.

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.