

**Anti-Complement factor B (human)  
 Mouse monoclonal antibody**

Subclass: IgG1/k

PRODUCT NO.

**HYB 008-06**

Clone: 6G11

PRESENTATION

Preparation: Protein-A/G purified

Content: Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.

Solvent: 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide

Storage: 4-8°C without exposure to light. No precautions necessary during handling.

ANTIGEN

Factor B is a component of the alternative pathway of complement activation. The protein consists of a single polypeptide chain of molecular mass 93 kDa, and has a concentration in plasma of approximately 200 µg/mL. Factor B is cleaved by factor D into the Ba component (30 kDa) and the Bb (63 kDa). Together with C3b, Bb forms the alternative pathway C3 convertase C3bBb (1).

IMMUNOGEN

Complement factor B purified from human plasma and adsorbed to aluminum hydroxide gel

SPECIFICITY

HYB 008-06 is specific for the Ba fragment of human factor B

EPI TOPE SPECIFICITY

Differs from that of HYB 008-02, HYB 008-04 and HAV 005-02.

REACTIVITY

HYB 008-06 reacts strongly with human factor B. A strong reaction is seen when coated directly onto the microtiter well and also when used as a biotinylated detection antibody in sandwich ELISA with a polyclonal capture antibody against factor B (e.g. Dako A0343). In Western blotting after SDS-PAGE, HYB 008-06 reacts with the Ba fragment of factor B in its non-reduced form only.

CULTURE MEDIUM

RPMI 1640 with 2-10% fetal calf serum

FUSION PARTNER

X63-Ag8.653

IMMUNIZATION

Female CF1 x BALB/c mice immunized by intraperitoneal injection

APPLICATION

Method	Usability	References
ELISA	Yes	
Immunoblotting	Yes	
Immunohistochemistry	Not determined	

REFERENCES

1. Kolb WP, Morrow PR, Tamerius JD (1989) Ba and Bb fragments of factor B activation: fragment production, biological activities, neopeptide expression and quantitation in clinical samples. Complement Inflamm 6:175-204.

**CONDITIONS**

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.