

**Anti-Glucagon-like peptide-1 (GLP-1, Mid-molecule specific)
 Mouse monoclonal antibody**

Subclass: IgG2a/k

PRODUCT NO.

HYB 147-08

Clone: 11E2

PRESENTATION

Preparation: Protein-A/G purified

Content: Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.

Solvent: 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide

Storage: 4-8°C without exposure to light. No precautions necessary during handling.

ANTIGEN

Glucagon-like peptide 1(7-36)amide (GLP-1(7-36)amide) is the principal active form of GLP-1, the other being GLP-1(7-37). GLP-1 is a peptide hormone of the glucagon family, produced by the L cells of the intestinal mucosa from the same prohormone as glucagon. The active forms are potent stimulators of glucose-dependent insulin secretion. The sequence of GLP-1 is fully conserved in all mammalian species examined so far.

IMMUNOGEN

Synthetic GLP-1(7-36)amide coupled to carrier and adsorbed onto aluminum hydroxide gel

SPECIFICITY

Reacts with all forms of GLP-1, including precursor and GLP-1(9-37) /GLP-1(9-36amide) metabolite

EPI TOPE SPECIFICITY

Mid-molecular epitope of GLP-1

REACTIVITY

HYB 147-08 binds to GLP-1 when coated directly onto the microtiter well, and binds GLP-1(7-36)amide in solution giving a K_a of 1.2×10^7 in inhibition ELISA. HYB 147-08 does not cross-react with coated glucagon or hGLP-2. In inhibition ELISA no binding of free glucagon in solution is detected, giving an estimated cross-reactivity of <0.2%. Although not tested, HYB 147-08 is likely to detect all known molecular forms of GLP-1 in immunohistochemistry.

In Western blotting a dilution guideline of 1/1000 has proved successful.

CULTURE MEDIUM

RPMI 1640 with 10% fetal calf serum

FUSION PARTNER

X63-Ag8.653

IMMUNIZATION

Female BALB/c mice immunized by intraperitoneal injections

APPLICATION

| Method | Usability | References |
|----------------------|----------------|------------|
| ELISA | Yes | |
| Immunoblotting | Yes | |
| Immunohistochemistry | Not determined | |

REFERENCES

.

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.