

**Anti-GLP-1(7-37) and GLP-1(7-36)amide (N-terminal specific)  
Mouse monoclonal antibody**

PRODUCT NO.

**ABS 033-04**

Subclass: IgG1/k

Clone: 4

PRESENTATION

Preparation: Protein-A purified

Content: Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.

Solvent: 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide

Storage: 4-8°C without exposure to light. No precautions necessary during handling.

ANTIGEN

Glucagon-like peptide 1(7-36)amide (GLP-1(7-36)amide) is the principal active form of GLP-1, the other being GLP-1(7-37). GLP-1 is a peptide hormone of the glucagon family, produced by the L cells of the intestinal mucosa from the same prohormone as glucagon. The active forms are potent stimulators of glucose-dependent insulin secretion. The sequence of GLP-1 is fully conserved in all mammalian species examined so far.

IMMUNOGEN

Synthetic GLP-1(7-17) adsorbed onto aluminum hydroxide gel

SPECIFICITY

ABS 033-04 binds the free N-terminus of GLP-1(7-37) and GLP-1(7-36)amide and shows <0.2% cross-reactivity with GLP-1(1-37), GLP-1(9-36)amide, glucagon, human GIP and exendin-4. Cross-reacts approximately 1% with human GLP-2.

EPI TOPE SPECIFICITY

Free N-terminus of GLP-1(7-37) and GLP-1(7-36)amide.

REACTIVITY

ABS 033-04 binds free GLP-1(7-37) and free GLP-1(7-36)amide in solution. The ELISA dilution guideline is to GLP-1(7-37) coated directly at 1 µg/mL. Not recommended for use in sandwich ELISA.

CULTURE MEDIUM

RPMI 1640 with 10% fetal calf serum

FUSION PARTNER

SP2mIL6

IMMUNIZATION

Female NMRixBALB/c mice immunized by intraperitoneal injection

APPLICATION

Method	Usability	References
ELISA	Yes (direct ELISA)	
Immunoblotting	Not determined	
Immunohistochemistry	Yes	

REFERENCES

**CONDITIONS**

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.